

Introduction to Metabolomics


10:00-11:00 Overview of techniques

Targeted and non-targeted metabolomics
(metabolite extraction procedures, equipment
GC-MS, HPLC-PDA-MS)

AIMS:


1. Experimental Design
2. Metabolite preparation: quenching and extraction
3. Metabolite separation and detection
4. Metabolite identification

The Metabolome



R. E. Gerszten & T. J. Wang (2008) Nature **451**, 949-952

1. Experimental design: things to consider before you start.



WHAT'S THE QUESTION?

Non-targeted metabolomics

comprehensive analysis of all the measurable analytes in a sample including chemical unknowns.

Hypothesis generating

Targeted metabolomics

the measurement of defined groups of chemically characterized and biochemically annotated metabolites

Hypothesis testing

1. Experimental design: technical considerations.

Material to be analysed:

- organism, tissue, organ, serum, cell, subcellular compartment?
- Contamination of sample with microbes
- How much sample 1-10mg MS
- 50-100mg NMR
- How many samples do you *need* for correct biological interpretation?
- How many samples do you *have* access to?



Sample preparation:

- How are you going to stop (**quench**) metabolism?
- Location – lab, field, hospital ward
- How are you going to extract the metabolites?



Analytical platforms:

What's available conduct the analysis?

- MS,
- HPLC,
- NMR
- GC



How are you going to **interpret** the data –

- univariate or
- multivariate statistics?




Will this answer your question?



AIMS:

1. Experimental Design
- 2. Metabolite preparation: quenching and extraction**
3. Metabolite separation and detection
4. Metabolite identification



Amino acids



O=C[C@H](O)[C@H](O)[C@@H](O)[C@H](O)CO
D-Glucose

O[C@H]1[C@H](O[C@H]2[C@@H](O)[C@H](O)[C@@H](O)[C@H]2O)[C@H](O)[C@@H](O)[C@H]1O

Raffinose

O=C(O)c1ccccc1O
Salicylic acid

O=C1OC(=O)C(O)C(O)CO1
Vitamine C

Fatty acids

Linoleic acid


Stearic acid

[illegible]

Lipids

Tricetyl glycerol

The diagram shows the chemical structure of Tricetyl glycerol, a triglyceride. It consists of a glycerol backbone esterified with three palmitic acid chains (16:0). The top chain is a long hydrocarbon tail. The middle chain is a long hydrocarbon tail with a carbonyl group (C=O) and an ester linkage (O-C=O). The bottom chain is a long hydrocarbon tail with a carbonyl group (C=O) and an ester linkage (O-C=O). The glycerol backbone is shown with a central carbon atom bonded to two oxygen atoms (O) and a hydrogen atom (H).



Solubility in water

2. Sample preparation: chemical and physical properties of metabolites

Molecular weight = the sum of weights of all atoms making the molecule,

$\text{H}_2\text{O} = 18$ (18 g per mol); lipids = >1000 g per mol

Molecular size = the 3D size of the structure, measured as Å

Polarity = differences in electronegativity:

Polar (large difference in positive and negative charges) (**hydrophilic**)

non(a)-polar compounds (no or little difference in charge) (**hydrophobic**)

more O and H = more polar; more N = less polar

Solubility = related to polarity, temperature and size


(like dissolves like)

To dissolve - particles need to separate and fit between the solvent spaces

eg, in polar metabolites, a positive end of a metabolite attaches to a negative end of solvent – cannot happen if a positive charge has no negative charge to attach to

Volatility = depends on boiling and melting point – liquid to gas phase
(more polar = less volatile)

Stability = thermal or oxidative instability



2. Sample preparation: Solvents for extraction

Polarity					
Highly non-polar (hydrophobic)			Highly polar (hydrophilic)		
Extract in Solvent	Heptane	Ethyl acetate		Acetone	Ethanol
	Hexane	Chloroform		Acetonitrile	Water
		Dichloromethane		Methanol	Perchloric acid
					NaCl
Metabolites extracted	Lipids	Carotenoids		Amino acids	Sugars
	Fatty acids	Chlorophylls		Organic acids	Nucleotides
	Waxes	Steroids		Organic amines	Phosphates
	Terpenes	Flavonoids		Alkaloids	

2. Sample preparation: Quenching

Quenching - Stop enzymatic metabolism

- Turnover rate is **fast**: reaction half lives < 1s
 - glucose to glucose-6-phosphate 0.3 to 1 mM per s
 - ATP used at a rate of 1.5mM per s
- Temperature:
 - Cold (< -40°C)
 - Hot (>80°C)
- pH
 - Acid (pH <2.0)
 - Alkaline (pH > 10)
- Common quenching methods:
 - Liquid nitrogen with or without freeze drying
 - Hot or cold Ethanol/Methanol
 - Perchloric acid
 - Sodium Hydroxide
 - Cold NaCl

Now metabolism is stopped, how do you extract metabolites?

2. Sample preparation: Extraction methods

Each group of metabolites will have an optimal extraction method (no single method will extract them all)

Liquid phase extraction

Grind sample, extract with solvent

Liquid : Liquid biphasic extraction

Take liquid extract, extract with another solvent

Solid : Liquid Extraction

Take liquid extract, extract with solid phase material

Volatile extraction

Steam distillation or headspace volatile collection

2. Sample preparation: Tissue disruption and Liquid phase extraction


Pestle and mortar



Grind tissue to fine powder in liquid nitrogen



Freeze dry material 1-2 days



Weigh a portion of freeze-dried sample and add extraction solvent



Sonicate for 20-30 minutes




Spin down tissue debris (requirement for cooled centrifuge?)



Recover supernatant for metabolite analysis



TissueLyser / Ball mill
Grinds frozen or freeze-dried material with or without extraction solvent



2. Sample preparation: Liquid: Liquid biphasic extraction

Biphasic partitioning allows polar and non-polar metabolites to be recovered separately

i) Extraction in polar solvent

extract in aqueous alcohols such as methanol, ethanol, and isopropanol.


Alcohols are miscible with water and may be used as aqueous phase co-solvents to enhance the solubility of less-polar metabolites during the extraction process

ii) Mix with non-polar solvent


organic solvents such as chloroform

CHCl_3 is volatile, nonreactive, immiscible with and denser than water

to form a two-phase system for liquid-liquid extraction



Sana and Fischer, Agilent Application Note,




2. Sample preparation: Common bi-phasic metabolite extraction procedure

Solvent mixture A = MeOH/CHCl₃/H₂O, 2.5:1:1, v/v/v at –20 °C;

Solvent mixture B = MeOH/CHCl₃, 1:1, v/v at –20 °C

Solvent C = deionised/distilled H₂O at 4 °C



What next – need to analyse metabolites

2. Sample preparation: Quality control

ALWAYS validate methodologies

- **Pool** from representative samples after extraction
- **Run** at the start and end, an every 5 or 10 samples during data acquisition
- Observe **technical** reproducibility
- **Spike** (add) extract with known amount of non-interfering substance – can you recover all of that spike after your analysis?

AIMS:

1. Experimental Design
2. Metabolite preparation: quenching and extraction
- 3. Metabolite separation and detection**
4. Metabolite identification

3. Metabolite separation


Separation Methods

- Thin layer chromatography
- High Performance Liquid Chromatography (HPLC)
- Gas chromatography (GC)
- Capillary electrophoresis (CE)


3. Metabolite separation : The basics

Thin Layer Chromatography using Paper or Silica Gel

- Aim is to separate (**resolve**) different metabolites in a mixture
- Maximum number of peaks that can be resolved is called '**peak capacity**'
- Can be increased by changing ratio of liquids/solvents or temperature




Black ink <http://www.teachengineering.org>



3. Metabolite separation : HPLC High Performance Liquid Chromatography

Same principle as TLC but **stationary phase** (resin) is packed tight in a narrow column and needs pumps to push the **mobile phase** (solvents) around the resin



Normal phase liquid chromatography





- column is packed full of a **polar** compound (eg. alkyl nitrile)
- non-polar mobile phase** such as hexane
 - good for lipids

Reverse-phase liquid chromatography

- column is packed with a **non-polar** silica compound
(eg C8 octasilane or C18 octadecylsilane)
- polar mobile** phases such as water/methanol/acetonitrile
- changes in pH, salts, solvent affect retention times
- good for phenolics, sugars, amino acids, drugs, pesticides

3. Metabolite separation : Reverse Phase HPLC


Different metabolites will have different retention times



Isocratic – same solvents ratios running through column eg. 100% Methanol

Gradient - change in solvent ratios over time
eg. start at 10 % MeOH
finish at 60% MeOH
over 20 minutes

3. Metabolite separation : HPLC The hardware



3. Metabolite separation : Detectors for HPLC

1. UV/VIS:

- Fixed wavelength
- Variable wavelength
- Diode array

2. Refractive index

3. Fluorescence

4. Conductivity

5. Antioxidant


6. Evaporative light scattering

7. Electrochemical

8. Mass Spectrometer

3. Metabolite separation : UV/VIS Detectors for HPLC

SPECTROPHOTOMETRY - Absorption of electromagnetic radiation



HPLC separation – set wavelength in UV range (284 nm)




Figure 1. Diagrammatic illustration of a UV-VIS detector optical system







Figure2. Diagrammatic illustration of a DAD optical system
Photodiode Array Detection




UV-absorbance Multi-scan wavelength (200 to 500nm)




What does this UV spectrum tell us about the compound?

Flavone – eg. Apigenin




Band II
240-280nm

Band I
at 304-350nm




Band I and II close

Flavonol – eg. Quercetin



Band II


Band I
at 352-385nm



Band I more defined and further away from band II

Characteristic UV absorption spectra can be used to i.d. flavonoids



Overlay unknown sample peak (eg, possible flavonol-glycoside)
with known standard peak




3. Metabolite separation : GC-FID

Gas Chromatography – Flame Ionisation Detection

Same principle as HPLC but use **GAS** rather than Liquid to separate metabolites



3. Metabolite separation : GC-FID



chemwiki.ucdavis.edu

Many polar compounds are NOT volatile
Sugars, amino acids, organic acids

Derivatisation – make compounds volatile
by making them more apolar


Silylation is the most widely used
technique.
acidic hydrogen replaced with with an
alkylsilyl group (eg. SiMe_3)
to form tri-methyl silyl (TMS) derivatives

Samples (about $1\mu\text{L}$) are injected into a
hot (250°C) glass tube where the
sample is vapourised

Vapour goes into the column


Separation based on the difference in
partition coefficients between the solid
(liquid) stationary phase and the
mobile gas phase

Increasing temperature biases
compounds to leave the stationary
phase and enter the gas phase



3. Metabolite separation : GC-FID

Detection: Flame Ionisation Detection



6. FID to give picogram-level detection

5. Ionization current between the two electrodes is directly proportional to the hydrocarbon concentration in the eluate that is burned by the flame


4. +ve ions migrate to the negative collector electrode and electrons migrate to the jet and are grounded

3. A high polarizing voltage is applied between the two electrodes

2. When sample emerges from column it is ionised in the flame producing +ve ions and electrons.

1. Mixture of Hydrogen and Air is used for flame. Burning hydrocarbon-free hydrogen in hydrocarbon-free air produces a negligible number of ions

Gas Chromatography (GC) – Flame Ionisation Detector (FID)



AIMS:

1. Experimental Design
2. Metabolite preparation: quenching and extraction
3. Metabolite separation and detection
- 4. Metabolite identification**

4. Metabolite identification

Mass Spectrometry

- Direct Injection Mass Spectrometry (DIMS)
- HPLC-PDA-MS
- GC-MS

4. Metabolite identification

What is Mass Spectroscopy?

Determination of the molecular weight of atoms and molecules

This is useful because the knowing the mass of a molecule is sometimes sufficient, often necessary and always useful in identifying it.

How is Mass Spectroscopy done?

A charge is put on the molecules of interest, (i.e., the analyte). Then measures how the trajectories of the resulting ions respond (in vacuum) to various combinations of electric and magnetic fields


What does Mass Spectroscopy measure?

measure the m/z ratio, not the mass


if an ion has multiple charges, the m/z will be significantly less than the actual mass

Knowing the mass of a molecule can help us to identify it.

Molecular mass or molecular weight refers to the mass of a molecule.
It is calculated as the sum of the mass of each constituent atom



4. Metabolite identification: The basic hardware



Ionisation

Individual metabolites in the sample are ionised and either become positively or negatively charged.

Acceleration (Separation)

These ions are then accelerated so that they all have the same amount of energy.

Deflection (Separation)

The ions are then deflected by a magnetic field according to their masses. The lighter and more charged they are, the more they are deflected.

Detection

The ions passing through the machine are detected electrically.

inlet

Ionisation
source

Mass
analyser

Ion
detector


4 basic components of all mass
spectrometers

4. Metabolite identification: Ionisation sources

The sine qua non of such a method (mass spectroscopy) is the conversion of neutral analyte molecules into ions.

For small and simple species the ionization is readily carried by gas-phase encounters between the neutral molecules and electrons, photons, or other ions.

*In recent years, the efforts of many investigators have led to **new techniques for producing ions of species too large and complex to be vaporized without substantial, even catastrophic, decomposition.***



http://en.wikipedia.org/wiki/John_Fenn


John B. Fenn, the originator of electrospray ionization for biomolecules and the 2002 Nobel Laureate in Chemistry
http://masspec.scripps.edu/mshistory/whatisms_details.php#Basics

4. Metabolite identification: Ionisation sources


Many different types:

- Electron Impact (EI) (hard)
- Atmospheric Pressure Chemical Ionisation (APCI) (good for polar compounds)
- Chemical Ionisation (CI)
- **Electrospray Ionisation (ESI) (good for polar compounds) (soft)**
- **Fast Atom Bombardment (FAB)**
- **Matrix Assisted Laser Desorption Ionisation (MALDI)**

ESI and MALDI are the ionisation sources of choice for biomolecular analysis.



4. Metabolite identification: Typical sensitivity and mass ranges allowed by different ionization techniques




http://masspec.scripps.edu/mshistory/whatisms_details


- electron ionization (EI),
- atmospheric pressure chemical ionization (APCI) while
- electrospray ionization (ESI), nanoelectrospray ionization (nanoESI),
- matrix-assisted laser desorption ionization (MALDI) have a high practical mass range.

4. Metabolite identification: Electrospray Ionisation (ESI)

evaporation of charged droplets



4. Metabolite identification: Electrospray Ionisation (ESI)



capillary

Sampling cone


4. Metabolite identification: Positive or Negative Ionisation

Positive ion mode: often by addition of H^+ , Na^+ , K^+

Negative ion mode: often by loss of H^+ or addition of Cl^-


If the sample has functional groups that readily accept a proton (H^+) then **positive** ion detection is used

e.g. amines $\text{R-NH}_2 + \text{H}^+ = \text{R-NH}_3^+$ as in proteins or peptides.



If the sample has functional groups that readily lose a proton then **negative** ion detection is used

e.g. carboxylic acids $\text{R-COOH} = \text{R-CO}_2^-$ and alcohols $\text{R-OH} = \text{R-O}^-$ as in saccharides or oligonucleotides



4. Metabolite identification: Separation



•The **mass analyser separates**, or **resolves**, the ions formed in the ionisation source of the MS according to their **mass-to-charge (m/z) ratios**.

•Quadrupole

- Four solid rods connected in parallel to a radio frequency (RF) generator and a **fixed** Direct Current (DC)
- Ions travel in between the four rods
- Trajectory of ion travel can be changed by changing RF.
- Unstable ions do not travel but selected stable ones do continue to travel to the detector.
- By scanning the RF field a broad m/z range (typically 100 to 4000) can be achieved in approximately one second.
- quadrupole analyzers typically scan up to m/z 3000

•Magnetic sectors –


- a charged molecule propelled through a magnetic field can be deflected by that field in a manner that depends on the mass and its charge



4. Metabolite identification: Separation Tandem Mass Analysis (MS/MS or MSⁿ)

Tandem mass analysis is used to sequence peptides, and structurally characterize carbohydrates, small oligo-nucleotides, and lipids

-A “triple-quad” is three of these put together in series.




4. Metabolite identification: MS/MS on Jasmonic acid

JA_100ng (2)

26/03/2020 17:14:00


JA 10 ng/ul

RT: 0.00 - 35.02 SM: 5B



NL: 2.46E7
m/z= 192.30-193.30 F:
+ c APCI corona Full
ms2 211.10@cid30.00
[100.00-220.00] MS
JA_100ng (2)

SIM
Selective ion
monitoring for
daughter fragment
ion from JA




$C_{12}H_{18}O_3$

Molecular Weight 210



JA_100ng (2) #2417 RT: 23.90 AV: 1 NL: 2.09E7

F: + c APCI corona Full ms2 211.10@cid30.00 [100.00-220.00]



← JA ion

4. Metabolite identification: Detection



The **detector** monitors the ion current, amplifies it and the signal is then transmitted to the data system where it is recorded in the form of **mass spectra** .



The **m/z values** of the ions are plotted against their **intensities** to show the **number of components** in the sample, the **molecular mass** of each component, and the **relative abundance** of the various components in the sample.

Accuracy: The ppm (parts per million) mass accuracy is a percent error quoting the difference between the measured and calculated mass for a particular ion. 0.1% would be equivalent to 1000 ppm error.

Standard 5ppm error is equivalent to 0.0005% - eg: for a mass 100 m/z the error would be +/- 0.05 and on mass 1000 error would be +/- 0.5

4. Metabolite identification: metabolic profiling


High Performance Liquid Chromatography (HPLC) – Photodiode array (PDA) – Mass spectrometry (MS)





Cyanidin 3-(3-malonyl glucoside)
Mr=534.90

448.9 (-malonyl)
287.04 (-hexose and malonyl)

High Performance Liquid Chromatography (HPLC) – Photodiode array (PDA) – Mass spectrometry (MS)



Metabolic Profiling - Gas Chromatography (GC) – Mass spectrometry (MS)





Unique fragmentation patterns

Matt Davey


GC-MS NIST fragment metabolite library – eg. fatty alcohol

Hit	SI	RSI	Prob	Name	Library Name
1	897	897	46.93	1-Octadecanol	MAINLIB
2	875	876	18.55	1-Nonadecene	MAINLIB
3	857	865	9.56	1-Pentadecanol	MAINLIB
4	850	856	7.32	1-Heptadecanol	MAINLIB
5	844	854	5.76	1-Tetradecanol	MAINLIB
6	838	886	4.52	1-Tridecanol	MAINLIB
7	818	832	2.06	1-Tetradecene	MAINLIB
8	814	834	1.74	1-Tridecene	MAINLIB
9	806	816	1.29	1-Pentadecene	MAINLIB
10	794	822	0.86	Cyclododecane	MAINLIB
11	782	817	0.57	1-Dodecanol	MAINLIB
12	766	840	0.33	1-Undecanol	MAINLIB
13	765	798	0.31	1-Dodecene	MAINLIB
14	734	777	0.08	1-Undecene	MAINLIB
15	710	758	0.02	1-Decanol	MAINLIB
16	690	790	0.00	1-Octadecanol	MAINLIB


1-Octadecanol
Formula C₁₈H₃₈O, MW 270, CAS# 112-92-5, Entry# 308
n-Octadecanol




Raw data - Library entry



Metabolite fingerprinting



Direct injection of crude extract



Could be mass/1 or mass/2

NMR
FTIR
RAMAN

**11:20-12:30 How to identify a
metabolite**
Online practical

1: How to identify metabolites from a mass spectrum

- KNApSack

Search by name of...

☒ Organism ☐ Metabolite

#> Enter name

List all

Search by molecular weight

#> Enter Mw ± margin

Mw ±

List

Search by molecular formula

#> Enter molecular formula

List

Molecular structure

Organism	Molecular formula	Metabolite	Mw
----------	-------------------	------------	----

Calculation mode for Mw ->>

Actual

Search by hierarchy

Search by Mass spectra

Display chart

Select MS data

Select MS data

Example

Arabidopsis1.txt

Arabidopsis2.txt


Arabidopsis3.txt

Arabidopsis4.txt

DIMS_example.txt

Latest update: 12/19 --2008 [Version KNApSack for Download -- v1.200.03]

...by Laboratory of Comparative Genomics, NAIST



KNApSack project was started on the April 1st in 2004 by 6 persons,
Ken Kurokawa, Yukiko Nakamura, Hiroko Asahi, Yoko Shinbo, Md. Altaf-Ul-Amin and Shigehiko Kanaya.
So we call this project KNApSack by picking and arranging the front letters of the family names of the members.
The pc means the Database can be accessed and used by a personal computer.

Current Status :

50897 Metabolites.
106418 Sp.-Metabolite pairs.

Statistics of genus

MassBank

<http://www.massbank.eu/MassBank/>

The screenshot shows the MassBank NORMAN MassBank Mass Spectral DataBase website. The browser window has a blue title bar with the text "MassBank | NORMAN MassBank Mass Spectral DataBase - Mozilla Firefox". The address bar shows the URL "http://www.massbank.eu/MassBank/". The website header features the MassBank logo and the text "High Resolution Mass Spectral Database". Below the header, the main content area is titled "NORMAN MassBank" and contains several search and browse options:

- Spectrum Search**: A search interface with a text input field and a "Search" button.
- Quick Search**: A search interface with a text input field and a "Search" button.
- Substructure Search**: A search interface with a chemical structure input field and a "Search" button.
- Browse Page**: A search interface with a text input field and a "Search" button.
- Peak Search**: A search interface with a text input field and a "Search" button.
- Spectral Browser**: A search interface with a text input field and a "Search" button.
- Record Index**: A search interface with a text input field and a "Search" button.

The website is displayed in a Mozilla Firefox browser window with multiple tabs open. The browser's address bar shows the URL "http://www.massbank.eu/MassBank/". The browser's status bar at the bottom shows the text "x".

•Kegg tab (www.genome.jp/kegg)

KEGG: Kyoto Encyclopedia of Genes and Genomes - Mozilla Firefox

File Edit View History Bookmarks Tools Help

KEGG: Kyoto Encyclopedia of Genes and Ge... +

www.genome.jp/kegg/

Property Details | Engl... Henderson-Hasselbalc... ChlamyCyc Phytosome v4.0: Sear... Chlamydomonas reinh... CyanoBase: Similarity ... Phytosome v4.0: Sear... Resources - GARNet TAIR - AraCyc



KEGG Search Help

> Japanese

KEGG Home
[Release notes](#)
[Current statistics](#)
[Plea from KEGG](#)

KEGG Database
[KEGG overview](#)
[Searching KEGG](#)
[KEGG mapping](#)
[Color codes](#)

KEGG Objects
[Pathway maps](#)
[Brite hierarchies](#)

KEGG Software
[KegTools](#)
[KEGG API](#)
[KGML](#)

KEGG FTP
[Subscription](#)

[GenomeNet](#)
[DBGET/LinkDB](#)
[Feedback](#)
[Kanehisa Labs](#)

KEGG: Kyoto Encyclopedia of Genes and Genomes

KEGG is a database resource for understanding high-level functions and utilities of the biological system, such as the cell, the organism and the ecosystem, from molecular-level information, especially large-scale molecular datasets generated by genome sequencing and other high-throughput experimental technologies (See [Release notes](#) for new and updated features).

- Main entry point to the KEGG web service**
 - KEGG2** [KEGG Table of Contents](#) [Update notes](#)
- Data-oriented entry points**
 - KEGG PATHWAY** [KEGG pathway maps](#) [[Pathway list](#)]
 - KEGG BRITE** [BRITE functional hierarchies](#) [[Brite list](#)]
 - KEGG MODULE** [KEGG modules](#) [[Module list](#)]
 - KEGG DISEASE** [Human diseases](#) [[Cancer](#) | [Infectious disease](#)]
 - KEGG DRUG** [Drugs](#) [[ATC drug classification](#)]
 - KEGG ORTHOLOGY** [Ortholog groups](#) [[KO system](#)]
 - KEGG GENOME** [Genomes](#) [[KEGG organisms](#)]
 - KEGG GENES** [Genes and proteins](#) [Release history](#)
 - KEGG LIGAND** [Chemical information](#) [[Reaction modules](#)]
- Entry point for wider society**
 - KEGG MEDICUS** [Health-related information resource](#)
- Organism-specific entry points**
 - KEGG Organisms** [Enter org code\(s\)](#) [Go](#) [hsa](#) [hsa eco](#)

Ionization method		Ionization source
Protonation	proton is added to a molecule, producing a net positive charge of 1+ for every proton added.	MALDI, ESI, and APCI
Deprotonation	Proton is removed from a molecule producing the net negative charge of 1-	MALDI, ESI, and APCI
Cationization	produces a charged complex by non-covalently adding a positively charged ion to a neutral molecule	MALDI, ESI, and APCI
Transfer of a charged molecule to the gas phase	transfer of compounds already charged in solution is normally achieved through the desorption or ejection of the charged species from the condensed phase into the gas phase	MALDI or ESI
electron ejection	ejection of an electron to produce a 1 ⁺ net positive charge	electron ionization (EI)
electron capture	a net negative charge of 1- is achieved with the absorption or capture of an electron	

